



QuickSleeper®

by **Dentalhitec**

More powerful, more ergonomic, more stylish,
opt for top-of-the-range digital anesthesia.

Intraosseous anesthesia, the standard for dental anesthesia

The secret of *QuickSleeper* lies in its ability to deposit the anesthetic solution as close as possible to the apices, where it is most effective.

By reducing the distance between the injection point and the target, you avoid excessive dilution of the anesthetic, guaranteeing **greater efficiency and flawless reliability**.



AN EFFICIENT PRACTICE !

Its effectiveness is based on a simple 3-step protocol :

1. Mucosal anesthesia with painless insertion.
2. Osseous perforation.
3. Injection close to the apices for immediate effectiveness.

Why opt for intraosseous anesthesia with *QuickSleeper*?

QuickSleeper offers the **most effective solution on the market for pain-free, stress-free procedures**. It is the only device that allows the anesthetic to be injected as close as possible to the apex, providing **immediate, profound and pain-free anesthesia for adults and children**.

Here is why dental surgeons love it :

- + Immediate, painless, profound anesthesia.
- + Effective even on mandibular molars with pulpitis.
- + No unnecessary numbing of soft tissue.
- + Reduced quantities injected.
- + Up to 6 teeth anesthetized with a single injection.
- + Enables treatments in multiple quadrants in the same visit.
- + Makes it easier to organize long appointments.



Discover the new *QuickSleeper6* !

QuickSleeper6 represents the culmination of three decades of continuous innovation, offering practitioners a tool that has been redesigned in terms of both ergonomics and performance.

This new version of the anesthetic process so familiar to practitioners represents a significant advance on previous generations, offering a connected device with enhanced ergonomics.

More than ever, *QuickSleeper6* improves the comfort of patients, practitioners and their teams in modern dental practices.

QUICKSLEEPER6 BOASTS A NUMBER OF KEY INNOVATIONS:

+ Optimized ergonomics : With a weight of just 185g and a length of 172mm, it is the world's most compact and lightweight electronic intraosseous anesthesia system.

+ Improved perforation capacity : Thanks to its unique electronically controlled needle rotation function, *QuickSleeper6* makes it easier to perforate the cortical bone, even in the most complex cases.

+ Customized injection : Mode + allows you to optimize the injection time for each patient.



More information on www.dentalhitec.com

The new *QuickSleeper6* is connected !



Total customization of injection speeds

Three injection curve options that can be adjusted from your smartphone or tablet. Do you prefer **patient comfort** or **saving time**? The choice is yours! Even for mucosal injections, **two default speeds** are available to suit every need.

Optimized ergonomics for unrivalled comfort

Manual or automatic delivery, set the parameters for your pedal with a single click. **Sounds can be customized**. Disable non-essential sounds or activate the "drop sonar" for perfect monitoring of mucosal injections, ideal for beginners with QuickSleeper6. **Quantity to be injected is programmable**, the device stops automatically at the desired quantity.

Analysis and recording of your practice

Detailed statistics : Consult your average quantities injected, the number of rotation cycles, and other parameters for precise monitoring of your procedures.

Activity history : detailed data (quantity, speed, rotation, etc.) is time-stamped for complete, transferable tracking of your anesthetics.

Direct access to Dentalhitec services : information, help, ordering, after-sales service.

All you need to know about Dentalhitec is just 1 click away.

Dentalhitec

Deliver the best anesthetics



ZI de l'Appentière · 49280 Mazières-en-Mauges · FRANCE
+33 (0)2 41 56 41 91 · mail@dentalhitec.com

www.dentalhitec.com

